

# CONTOUR LEGACY Leg Pillow




## CUSHION YOUR UPPER THIGHS AND KNEES

Longer support to keep your hips, legs, and knees aligned – extra thick support at the top to cushion your knees - no more rubbing bone on bone.



## MOVE FREELY, SLEEP IN ANY POSITION

Tapered edges and concave center cradles your thighs and knees snugly allowing you to move freely, change sleep positions, and bend your knees.



## STAY COMFORTABLE, STAY COOL

Super soft memory foam keeps its shape and is ventilated to stay cool.



## BREATHABLE, WASHABLE COVER

Plush, soft cover easily removes for washing.



## ENJOY FULL SUPPORT HIP TO ANKLES

With our specially made case designed to hold TWO Legacy Leg Pillows together.

The Legacy Pillow will give you the best night's sleep guaranteed. The brand-new shape hugs your thighs, stays in place, and provides more comfort than ever before. The Legacy pillow is extra-long for more support and has ventilated holes to help keep you cool at night. It will not crowd your pelvic area due to the crescent cutout at the top and the extra thick foam on the bottom ensures legs are kept in alignment and do not drop forward.

### Legacy Pillow Use:

Place Legacy pillow between legs with thicker portion on knees. The Legacy pillow should feel natural between the legs and hug thighs comfortably. Cleaning: Remove your Legacy pillow cover by unzipping it. Remove the cover from the foam and wash the cover on a gentle cycle and then dry it on low heat.



### Double Legacy Pillow Use:

Place first Legacy pillow between thighs with thicker cushioning on knees as normal. Then place the second pillow between your calves with the thicker cushion facing the knees. The second pillow can be adjusted to a position that best suits your comfort.

We are so confident you will love your purchase that we now offer a full 90 Day Comfort Guarantee on the Legacy Leg Pillow and will waive all restocking and disposal fees.




# CONTOUR LEGACY Leg Pillow

---

## **Optional Legacy Pillow Case:**

Slightly compress legacy pillow to fit through case opening with the tapered side being inserted first. Then adjust case as needed so that the case seams line up with pillow edges. For single pillow use, snap the end fastener together to keep it tucked away while sleeping.



If using two Legacy pillows and cases, snap one case to the other to create a full-length leg support case.



**Cleaning:** Remove the case from the pillow and place it in the washing machine on a gentle cycle and then dry it on low heat.

We are so confident you will love your purchase that we now offer a full 90 Day Comfort Guarantee on the Legacy Leg Pillow Case and will waive all restocking and disposal fees.

Contour is a registered trademark under license from Banyan Licensing, L.L.C © 2018. All rights reserved.

---

## **How to reduce the odor of the foam?**

This product is made from all new materials and, like a new car, will have a new product smell. You may reduce the smell by removing your Legacy Leg Pillow from the package to air out. Sometimes this can take a few days. To completely remove the “new smell,” the cover may also be laundered. Please follow instructions on the care tag.

## **What causes foam “yellowing”?**

Exposure to fresh air and light may cause the foam to “yellow.” This is normal and is more visible on white foam.